

TEST STRIPS FOR PERACETIC ACID (PAA)

These test strips are the widest range available on the market. They are designed for ranges of PERASAN® peracetic acid from 10-160 ppm. They are convenient and come packaged in an air-tight plastic vial of 50 test strips each. Simply dip the test strip and read the appropriate intensity of the color strip within 15-30 seconds. Readings taken beyond 30 seconds are unreliable.


These strips can be used quite effectively, especially in systems that are very high in salt brine and other electrolytes. In this case most test procedures involving thiosulfate-based titration are quite unreliable and inaccurate and will lead to false incorrect high concentrations. This is true even for Enviro Tech's titration procedure, which is otherwise regarded as the most reliable on the market.

Most organic matter, electrolyte and salt/brine solutions do not affect these test strips or their accuracy, so they can be relied upon to give reasonably accurate and quick results for the majority of use applications. Remember that test strips of this type will NOT give complete color matches. One must depend on the relative depth of shading by comparing the two closest matches, and the actual result will undoubtedly fall within the "range" between the two color shades, which is usually accurate enough for evaluation purposes. More precise concentrations can be achieved using Enviro Tech's drop test procedure, except as noted in the preceding paragraph.

PACKAGING: Vials of 50 strips. Available in individual vials or sets of 10 vials.

SHIPPING: From Modesto, CA, USA

ORDERING: 209-581-9576 or visit us @ envirotechinc.com



PERASAN®
0-160ppm
50 Test Strips
Peracetic Acid

Enviro Tech

Peracetic Acid ppm

0
10
20
40
60
85
160

TO USE:
Using a cup-size sample, immerse test strip and remove with pads face up. **DO NOT SHAKE OFF EXCESS WATER.** Wait 10 seconds and compare to color chart immediately. Discard strip after 30 seconds.

TIPS

- Keep wet fingers out of the vial.
- Close vial tightly after removing strip.
- Store in a cool, dry place.

USA

Enviro Tech
Modesto, CA 95358
www.envirotechinc.com
Ph. 209-581-9576

Code 3000ET Lot 4410 Exp. 6/06